



## TACTICAL MECHANICAL TOURNIQUET

### DESIGNED BY THE US MILITARY

Extremity hemorrhage is a major cause of preventable death globally and the number-one cause of preventable death on the battlefield.<sup>1</sup> Death from blood loss can occur in as little as five minutes from a severe extremity injury. Tourniquets are a proven device to stop severe, life-threatening bleeding.

Developed by the US military and recommended by the Committee on Tactical Combat Casualty Care (CoTCCC), the TMT™ is used by militaries, medical professionals, and public safety agencies worldwide.

1. Eastridge, Brian J et al. "Death on the battlefield (2001-2011): Implications for the Future of Combat Casualty Care." *The Journal of Trauma and Acute Care Surgery* vol. 73,6 Suppl 5 (2012): S431-7. doi:10.1097/TA.0b013e3182755dcc

**NO INTERNAL BAND:** Eases slack removal and allows for rapid resetting.

**2" WIDE:** Achieves occlusion at lower pressures.

**UV-PROTECTIVE ADDITIVE TO POLYMERS:** Durable enough to be stored outside or packaging.

**TWO LOCKING MECHANISMS:** Double-routed speed buckle and hook and loop backed strap prevents slippage and pressure loss.

**QUICK-DISCONNECT BUCKLE:** Prevents misrouting of tail and allows for trapped-limb application.

**ENHANCED TORSION BAR SECURING:** Audible signature "click" identifies the torsion bar is secure.

**MULTIFUNCTION:** Undeniable advantage when combined with the TQ-RAM† for point pressure over wound packing and additionally provides proximal hemorrhage control.

#### SPECIFICATIONS

<b>PN</b>	PDF 185 (TMT Black [UK/US]) PDF 187 (TMT International)	31-115 (Berry Compliant)
<b>NSN</b>	6515-99-258-6320	6515-01-656-6191
<b>CE marked</b>	✓	✗
<b>Purchase</b>	Direct	DAPA, GSA, VA DAPA
<b>Packaged dimensions</b>	4.5" x 2.5" x 2" (11.4 x 6.4 x 5.1cm)	4.5" x 2.5" x 2" (11.4 x 6.4 x 5.1cm)
<b>Weight</b>	3.5oz (99.2g)	3.5oz (99.2g)
<b>Shelf life*</b>	7 years	10 years
<b>Misc.</b>	Accommodates limb circumferences 5.75"-38" (14.6cm-96.5cm)	Accommodates limb circumferences 5.75"-38" (14.6cm-96.5cm)

\* From date of manufacture

† May not be available in your region.



visit [safeguardmedical.com](http://safeguardmedical.com)  
 call (US) +1 (855) 428-6074 / (Int'l) +44 (0) 1568 613942  
 email (US) [info@safeguardmedical.com](mailto:info@safeguardmedical.com) / (Int'l) [sales.admin@safeguardmedical.com](mailto:sales.admin@safeguardmedical.com)